







Power Wheelchair Mini Proportional Joystick Comparison Chart

Joystick Name	Photo	Manufacturer	Force	Travel	Isometric	Switch Ports	Electronic Packages	Programmable	Sealed/ Moisture resistant
Molecule		Adaptive Switch Labs (ASL)	8 grams	3 mm	No	On in-line box	LiNX MK6i Q-Logic R-net	Yes	No
Micro Extremity Control (MEC)		ASL	18 grams	13 mm	No	Push down on handle for Reset	Mk6i Q-Logic R-net	Can make Reverse act as Reset	No
Extremity Control Joystick		ASL	120 grams	6 mm	No	None	Mk6i Q-Logic R-net	Can make Reverse act as Reset	Yes
Micro Mini		ASL	Not measured by ASL	minimal	Yes	None	Mk6i Q-Logic R-net	No	No
Mini Joystick		HMC/ Permobil	13 grams	3 mm	No	1 jack on interface box for Reset	Mk6i Q-Logic R-net	No	No
Compact Joystick - Lite		HMC/ Permobil	50 grams	17 mm	No	2 jacks on joystick, Power and Reset	Mk6i Q-Logic R-net	No	No

Joystick Name	Photo	Manufacturer	Force	Travel	Isometric	Switch Ports	Electronic Packages	Programmable	Sealed/ Moisture resistant
Multi Joystick		Stealth/ Mo-Vis	50 grams	8 mm	No	2 jacks on joystick and interface box, Power and Reset	Q-Logic R-net	Yes	No
Micro joystick		Stealth/ Mo-Vis	8.5 grams	3.3 mm	No	2 jacks on interface box, Power and Reset	Q-Logic R-net	Yes	No
Precision Mini Proportional (PMP) joystick		Stealth iDrive	43 grams	7 mm	No	On iDrive CPU (optional)	Mk6i Q-Logic R-net	Yes (also non-programmable version)	Yes
MicroPilot		Switch-it/ Sunrise	11-50 grams adjustable	0 mm	Yes	None	Mk6i Q-Logic R-net	No	No
MicroGuide		Switch-it/ Sunrise	25 grams	6 mm	No	None	Mk6i Q-Logic R-net	No	No

Comparison information: Standard Joystick Force 250 grams, Standard Joystick Throw 28 mm